



1220 E. Birch St. #201 Brea, CA 92821 (800)750-5004

Haupt Dental Lab Inc.

## Denture Conversion Materials List

### Equipment

- Handheld Light Unit
- Bench Lathe
- Electric Lab Handpiece
- Dust Hood with Vacuum
- Plastic Alcohol Torch (not butane torch) - Keystone

### Hand Tools

- Rubber Dam Punch
- Bard Parker Scalpel Handle
- #25 Redwood Plaster Knife #6/7
- Small Surgical Scissors
- Perio Probe / Explorer
- Prosthetic Unigrip Screw Drivers - Nobel
- Prosthetic Torque Driver - Nobel

### Supplies

- Carbide Burs
  - H257 AEF - Komet
  - H71E - (014) - Komet
- Lathe Brush Wheels (Chung King Set - Keystone)
- Lathe Pumice/Rag Wheels - Keystone
- Medium Grit Pumice - Any
- Acrylic Polish Compound - "Shure Shine Plate Polish" - Keystone
- Implant/Abutment Replicas - Nobel Biocare
- Implant/Abutment Protection Caps - Nobel Biocare
- 20mm Guide Pins - Nobel Biocare
- PVS Rigid Bite Registration - Any
- PVS Lab Putty 80-90 Shore Hardness - Any
- Fast Set Acrylic (6 minute) - Lang
- Light Cure Resin (pink) - Triad
- Self Cure Composite - "Quick Up" VOCO