

Denture



Full Denture

Our full dentures provide an economical solution for full-arch replacement. Utilizing hand-made, high impact denture teeth and pressure moulded acrylics, our dentures look great and fit extremely well. Additional reinforcements can also be designed to integrate within the acrylic for superior strength and durability.



Palatal View



LOCATOR®

Denture



Locator Denture

Dental implants provide increased stability over traditional dentures. Using the Patented Locator attachment, patients can have up to 5lbs. of vertical resiliency per/implant on their denture. This increase in stability, affords the elimination of the denture palate on most cases.

These dentures look and fit extremely well. Additional reinforcement can also be designed to integrate within the acrylic for superior strength and durability.



2-4 Implants



Palatal View



ALL ON 4[©]

Hybrid



All on 4 Hybrid

Maximum function and minimal surgery! The All on 4 treatment concept has been proven to provide the best of both. The all on 4 hybrid denture sometimes referred to as "teeth in a day," provides the patient with immediate function in a fixed bridge solution.

Combining custom hand-made resin teeth and CAD/CAM titanium bars, the All on 4 hybrid final bridge is a masterpiece of dental material collaboration.



4 Implants



CAD/CAM Bar



FIXED ZIRCONIA PORCELAIN

Bridge



Fixed Zirconia/e.Max Bridge

The finest materials, hand built multi-layered porcelain, unsurpassed esthetics. Implant supported implant fixed porcelain bridges are fabricated using the latest in CAD/CAM zirconia technology and custom milled high noble alloys. Dental porcelains are then applied with the mastery of over 40 years of artistic talents to create the most life-like restorations.



4-6 Implants



Zirconia Framework



FIXED PORCELAIN

Bridge



Fixed Porcelain Bridge

The finest materials, hand built multi-layered porcelain, unsurpassed esthetics. Implant supported implant fixed porcelain bridges are fabricated using the latest in CAD/CAM zirconia technology and custom milled high noble alloys. Dental porcelains are then applied with the mastery of over 40 years of artistic talents to create the most life-like restorations.



6-8 Implants



CAD/CAM Abutments



MILLED BAR

Denture



Milled Bar Denture

Combining a "fixed" bridge with a removable denture affords options. Many patients with existing implants are looking for maximum stability and retention of their removable denture. With CAD/CAM precision milling technology, a titanium bar can be fabricated to fit the existing dental implants. The overdenture is then created to restore the missing teeth and soft tissues.

Within the overdenture, a precision male/female locking attachment is integrated to lock the prosthesis in place offering absolute rigidity.



4-6 Implants



CAD/CAM Bar

